



Lullaby easy handle  
Medical handles

C O N E A S Y

**MOQ:** 3.000 pcs  
**Colours:** Transparent  
**Material:** PC - Medical Polycarbonate  
**Description**  
Lullaby easy handle with lateral wings “easy”

**Note**  
Other colours and material available on request

| COD.      | GAUGES             | ØCANULA | COLOURS     | NOTE                      |
|-----------|--------------------|---------|-------------|---------------------------|
| CONEASY14 | 14G                | 2,00 mm | Transparent | Doesn't need plastic bush |
| CONEASY15 | 15G                | 1,80 mm | Transparent | Doesn't need plastic bush |
| CONEASY16 | 16G                | 1,60 mm | Transparent | Doesn't need plastic bush |
| CONEASYCS | FROM 17G<br>TO 27G |         | Transparent | *see plastic bush         |