



Cap with thread

T A P F I

MOQ: 3.000 pcs

Colours: Different colours

Material: LDPE – low density polyethylene

Description

Cap with thread

Note

Other colours and material available on request

COD.	GAUGES	ØCANULA	COLOURS	NOTE
TAPFI14G	14G	2,00 mm	Dark green	
TAPFI14CS	14GC	2,10 mm	Dark green	
TAPFI16G	16G	1,60 mm	White	
TAPFI18G	18G	1,27 mm	Pink	
TAPFI18CS	18GC	1,52 mm	Pink	
TAPFI20G	20G	0,90 mm	Yellow	
TAPFI20C	20GC	1,15 mm	Yellow	