



## Anatomic handle without luer Medical handles

CONANASL

**MOQ:** 3.000 pcs

**Colours:** Transparent

**Material:** PC - Medical Polycarbonate



### Description

Anatomic handle for biopsy needles

### Note

Other colours and material available on request


COD.	GAUGES	ØCANULA	COLOURS	NOTE
CONANASL	0,75		Transparent Blue	



MATERIAL	CODE
PC	CONANACSxx
MAKROLON	CONANACSxx M
ZYLAR	CONANACSxx Y

CONANASL	WITHOUT LUER
CONANACSS	WITH THREAD
CONANACSF	WITHOUT THREAD
CONANACS	UNIVERSAL
CODE	DESCRIPTION


Changes	Signature	Date
3		
2		
1		

 **CARRO LUIGI Srl**  
 25030 RUDIANO (BS)  
 Via Lavoro & Industria

**M0803.04 rev.00**

General Tolerances **ISO 2768-mk**

0.5 < x < 6	+0.1
6 < x < 30	+0.2
30 < x < 120	+0.3
120 < x < 400	+0.5

<u>Mould Denomination</u>		<u>N° Cavities</u>		<u>Order</u>	
UNIVERSAL ANATOMIC HUB		---		---	
<u>Member Denomination</u>		<u>Costumer code</u>		---	
UNIVERSAL ANATOMIC HUB		---		---	
<u>Date</u>	<u>Molding Material</u>	<u>C.tion</u>	<u>Machine</u>	<u>Mould Weight</u>	<u>Costumer part</u>
02/12/2022	SEE TABLE	---	---	---	SEE TABLE
<u>Drawn By</u>	<u>Approved By</u>		<u>Eyebolts</u>	<u>Mould Dimensions</u>	<b>A4</b> Sheet
Cassani M.	Direzione C.A.	---	---	---	
<u>N°Par.</u>	<u>N°Pez.</u>	<u>Dimensions</u>	<u>Material</u>	<u>Treatment</u>	<u>Scale</u>
---	---	---	---	---	-/-