



## Octagonal handle

### Medical handles

C A N A O T

**MOQ:** 3.000 pcs

**Colours:** Transparent

**Material:** PC - Medical polycarbonate

**Description**

Octagonal anatomic handle for biopsy needles

**Note**

Other colours and material available on request

COD.	GAUGES	ØCANULA	COLOURS	NOTE
CANAOT14	14G	2,00 mm	Transparent	Doesn't need plastic bush
CANAOT15	15G	1,80 mm	Transparent	Doesn't need plastic bush
CANAOT16	16G	1,60 mm	Transparent	Doesn't need plastic bush
CANAOTCS	FROM 17G TO 27G		Transparent	*see plastic bush