



## Anatomic handle

### Medical handles

C O N A N A

**MOQ:** 3.000 pcs

**Colours:** Transparent

**Material:** PC - Medical Polycarbonate

**Description**

Anatomic handle for biopsy needles

**Note**

Other colours and material available on request

COD.	GAUGES	ØCANULA	COLOURS	NOTE
CONANA14	14G	2,00 mm	Transparent	Doesn't need plastic bush
CONANA15	15G	1,80 mm	Transparent	Doesn't need plastic bush
CONANA16	16G	1,60 mm	Transparent	Doesn't need plastic bush
CONANACS	FROM 17G TO 27G		Transparent	*see plastic bush