



Lullaby easy handle Medical handles

C O N E A S Y

MOQ: 3.000 pcs

Colours: Transparent

Material: PC - Medical Polycarbonate

Description

Lullaby easy handle with lateral wings "easy"

Note

Other colours and material available on request

COD.	GAUGES	ØCANULA	COLOURS	NOTE
CONEASY14	14G	2,00 mm	Transparent	Doesn't need plastic bush
CONEASY15	15G	1,80 mm	Transparent	Doesn't need plastic bush
CONEASY16	16G	1,60 mm	Transparent	Doesn't need plastic bush
CONEASYCS	FROM 17G TO 27G		Transparent	*see plastic bush