



Lullaby revolving handle

Medical handles

CONREV

MOQ: 3.000 pcs

Colours: Transparent

Material: PC - Medical Polycarbonate

Description

Lullaby revolving handle for biopsy needles with wings and revolving ferrule

Note

Other colours and material available on request

COD.	GAUGES	ØCANULA	COLOURS	NOTE
CONREV14	14G	2,00 mm	Transparent	Doesn't need plastic bush
CONREV15	15G	1,80 mm	Transparent	Doesn't need plastic bush
CONREV16	16G	1,60 mm	Transparent	Doesn't need plastic bush
CONREVC	FROM 17G TO 27G		Transparent	*see plastic bush