


YELLOW	NGB20G	20G	Ø0.90
PINK	NGB18G	18G	Ø1.30
WHITE	NGB16G	16G	Ø1.70
LIGHT GREEN	NGB14G	14G	Ø2.20
RED	NGB13G	13G	Ø2.45
DARK GREEN	NGB11G	11G	Ø2.95
COLOR	CODE	GAUGE	ØA

3	2	1
Signature		
Date		

Changes



CARRO LUIGI Srl
 25030 RUDIANO (BS)
 Via Lavoro & Industria

M0803.04 rev.00

General Tolerances **ISO 2768-mk**

0.5 < x < 6	+0.1
6 < x < 30	+0.2
30 < x < 120	+0.3
120 < x < 400	+0.5

Mould Denomination

N° Cavities

Order

Member Denomination

Customer code

NEEDLE GUIDE

Date
17/01/2023

Molding Material
ABS

C.tion

Machine

Mould Weight

Customer part

Drawn By
Cassani M.

Approved By
Direzione C.A.



Eyebolts

Mould Dimensions

SEE TABLE

N°Par. N°Pez.

Dimensions

Material

Treatment

Scale

A4

Sheet

1/1