



## Plastic bush

### Medical connectors

P B \_ \_ G

**MOQ:** 3.000 pcs

**Colours:** Transparent

**Material:** Acrylic

#### Description

Plastic bush to define the internal diameter of the universal handles

#### Note

To combine with universal handles

COD.	GAUGES	ØEXTERNAL CANULA	COLOURS	NOTE
PB17G	17G	1,47 mm	Transparent	To insert in universal handles (all models)*
PB18G	18G	1,27 mm	Transparent	*
PB19G	19G	1,10 mm	Transparent	*
PB20G	20G	0,90 mm	Transparent	*
PB21G	21G	0,80 mm	Transparent	*
PB22G	22G	0,70 mm	Transparent	*
PB23G	23G	0,60 mm	Transparent	*
PB24G	24G	0,55 mm	Transparent	*
PB25G	25G	0,50 mm	Transparent	*
PB26G	26G	0,45 mm	Transparent	*
PB27G	27G	0,40 mm	Transparent	*