



Mini Stopper

Medical stoppers

S T 2 1 1 _

MOQ: 3,000 pcs

Colours: Different colours

Material: Medical ABS

Description: Sliding stoppers for biopsy needles, used as indicator of penetration depth of the needle.

Note: Other diameters available on request

| COD. | GAUGES | ØCANULA | COLOURS | NOTE |
|----------|--------|---------|-------------|------|
| ST211_14 | 14G | 2,00 mm | Dark green | |
| ST211_15 | 15G | 1,80 mm | Blue | |
| ST211_16 | 16G | 1,60 mm | White | |
| ST211_17 | 17G | 1,46 mm | Red | |
| ST211_18 | 18G | 1,26 mm | Pink | |
| ST211_19 | 19G | 1,10 mm | Beige | |
| ST211_20 | 20G | 0,90 mm | Yellow | |
| ST211_21 | 21G | 0,80 mm | Light green | |
| ST211_22 | 22G | 0,70 mm | Black | |