



Straight connector tube-thread

Medical connectors

S T 2 3 0

MOQ: 3.000

Colours: Bianco

Material: Nylon

Description

Straight connector tube-thread

Note

Other colours and material available on request

COD.	GAUGES	ØCANULA	COLOURS	NOTE
230-3214	3,2 mm	/	White	Thread 1/4
230-4014	4,0 mm	/	White	Thread 1/4
230-4814	4,8 mm	/	White	Thread 1/4
230-6414	6,4 mm	/	White	Thread 1/4
230-2418	2,4 mm	/	White	Thread 1/8
230-3218	3,2 mm	/	White	Thread 1/8
230-4018	4,0 mm	/	White	Thread 1/8
230-4818	4,8 mm	/	White	Thread 1/8