



Straight connector

Medical connectors

S T 2 4 0/2 5 0

MOQ: 3.000 pcs

Colours: White

Material: Nylon

Description

Straight connector

Note

Other colours and material available on request

COD.	GAUGES	ØCANULA	COLOURS	NOTE
240_4	4 mm	/	White	
240_48	4,8 mm	/	White	
240_64	6,4 mm	/	White	
250_16	1,6 mm	/	White	
250_24	2,4 mm	/	White	
250_32	3,2 mm	/	White	