



## Stopper

### Medical stoppers

**S T O P P**

**MOQ:** 3.000 pcs

**Colours:** Different colours

**Material:** Medical ABS

**Description:** Sliding stoppers for biopsy needles, used as indicator of penetration depth of the needle.

**Note:** Other diameters available on request

COD.	GAUGES	ØCANULA	COLOURS	NOTE
STOPP14G	14G	2,00 mm	Dark green	
STOPP15G	15G	1,80 mm	Blue	
STOPP16G	16G	1,60 mm	White	Coaxial measure available
STOPP17G	17G	1,47 mm	Red	
STOPP18G	18G	1,27 mm	Pink	Coaxial measure available
STOPP19G	19G	1,10 mm	Beige	Coaxial measure available
STOPP20G	20G	0,90 mm	Yellow	Coaxial measure available
STOPP21G	21G	0,80 mm	Light green	
STOPP22G	22G	0,70 mm	Black	Coaxial measure available
STOPP23G	23G	0,60 mm	Light blue	
STOPP24G	24G	0,55 mm	Violet (Mauve)	
STOPP25G	25G	0,50 mm	Orange	
STOPP26G	26G	0,45 mm	Brown	
STOPP27G	27G	0,40 mm	Grey	