



## Rectangular cap Medical caps

T A P R E T

**MOQ:** 3.000 pcs

**Colours:** Different colours

**Material:** Medical ABS

**Description:** Rectangular cap, suitable for rectangular and anatomic handles.

Related: Rectangular and anatomic handles, stoppers and wings.

COD.	GAUGES	ØCANULA	COLOURS	NOTE	
TAPRET14	14G	1,70 mm	Dark green	Rectangular and anatomic handles	Coaxial measure available
TAPRET15	15G	1,50 mm	Blue	Rectangular and anatomic handles	
TAPRET16	16G	1,30 mm	White	Rectangular and anatomic handles	Coaxial measure available
TAPRET17	17G	1,10 mm	Red	Rectangular and anatomic handles	
TAPRET18	18G	1,00 mm	Pink	Rectangular and anatomic handles	Coaxial measure available
TAPRET19	19G	0,80 mm	Beige	Rectangular and anatomic handles	Coaxial measure available
TAPRET20	20G	0,65 mm	Yellow	Rectangular and anatomic handles	Coaxial measure available
TAPRET21	21G	0,55 mm	Light green	Rectangular and anatomic handles	
TAPRET22	22G	0,45 mm	Black	Rectangular and anatomic handles	
TAPRET23	23G	0,35 mm	Light blue	Rectangular and anatomic handles	
TAPRET24	24G	0,35 mm	Violet (Mauve)	Rectangular and anatomic handles	
TAPRET25	25G	0,30 mm	Orange	Rectangular and anatomic handles	
TAPRET26	26G	0,30 mm	Brown	Rectangular and anatomic handles	
TAPRET27	27G	0,25 mm	Grey	Rectangular and anatomic handles	