



PC	//
MAKROLON	CONFEMxx+M

CONFEM	GAUGE	ØA	ØB	ØC	COLOUR	NOTE	
CONFEM22	22G	0.70	0.40	1.10	BLACK	NOT COVERING BLACK	
CONFEM21	21G	0.85	0.55	1.15	LIGHT GREEN	NOT COVERING L. GREEN	
CONFEM20	20G	0.93	0.60	1.25	YELLOW	NOT COVERING YELLOW	
CONFEM19	19G	1.10	0.80	1.50	BEIGE	NOT COVERING BEIGE	
CONFEM18	18G	1.32	1.10	1.60	PINK	NOT COVERING PINK	
CONFEM17	17G	1.50	1.14	1.99	RED	NOT COVERING RED	
CONFEM16	16G	1.63	1.30	1.99	WHITE	NOT COVERING WHITE	
CONFEM16M					WHITE	NOT COVERING WHITE	
CONFEM15	15G	1.83	1.54	1.99	TRANSPARENT		
CONFEM14	14G	2.04	1.75	2.04	TRANSPARENT		
	CODE	GAUGE	ØA	ØB	ØC	COLOUR	NOTE

3	2	1
---	---	---

CARRO LUIGI Srl
 25030 RUDIANO (BS)
 Via Lavoro & Industria

M0803.04 rev.00

General Tolerances ISO 2768-mk

0.5 <x<6	+0.1
6 <x<30	+0.2
30 <x<120	+0.3
120 <x<400	+0.5

Mould Denomination _____				N°Cavities ---		Order ---				
Member Denomination FEMALE LUER HANDLE				Customer code ---						
Date 09/11/2022		Molding Material SEE TABLE		Machine _____		Customer part SEE TABLE				
Drawn By Cassani M.		Approved By Direzione C.A.		Eyebolts _____		Mould Dimensions _____				
N°Par. _____		N°Pez. _____		Dimensions _____		Material _____				
Treatment _____		Scale _____		<table border="1"> <tr> <td rowspan="2">A4</td> <td>Sheet</td> </tr> <tr> <td>1/1</td> </tr> </table>				A4	Sheet	1/1
A4	Sheet									
	1/1									

Changes

Signature

Date