



PC	//
MAKROLON	CONFEMxx+M

CONFEM	GAUGE	ØA	ØB	ØC	COLOUR	NOTE
CONFEM22	22G	0.70	0.40	1.10	BLACK	NOT COVERING BLACK
CONFEM21	21G	0.85	0.55	1.15	LIGHT GREEN	NOT COVERING L. GREEN
CONFEM20	20G	0.93	0.60	1.25	YELLOW	NOT COVERING YELLOW
CONFEM19	19G	1.10	0.80	1.50	BEIGE	NOT COVERING BEIGE
CONFEM18	18G	1.32	1.10	1.60	PINK	NOT COVERING PINK
CONFEM17	17G	1.50	1.14	1.99	RED	NOT COVERING RED
CONFEM16	16G	1.63	1.30	1.99	WHITE	NOT COVERING WHITE
CONFEM16M					WHITE	NOT COVERING WHITE
CONFEM15	15G	1.83	1.54	1.99	TRANSPARENT	
CONFEM14	14G	2.04	1.75	2.04	TRANSPARENT	
CODE	GAUGE	ØA	ØB	ØC	COLOUR	NOTE

3	2	1
---	---	---

CARRO LUIGI Srl
 25030 RUDIANO (BS)
 Via Lavoro & Industria

M0803.04 rev.00

General Tolerances **ISO 2768-mk**

0.5 < x < 6	+0.1
6 < x < 30	+0.2
30 < x < 120	+0.3
120 < x < 400	+0.5

<u>Mould Denomination</u>					<u>N° Cavities</u>		<u>Order</u>	
<u>Member Denomination</u>					<u>Costumer code</u>		<u>Costumer part</u>	
FEMALE LUER HANDLE					<u>SEE TABLE</u>		<u>SEE TABLE</u>	
<u>Date</u>	<u>Molding Material</u>	<u>C.tion</u>	<u>Machine</u>	<u>Mould Weight</u>	<u>Costumer part</u>			
09/11/2022	SEE TABLE				<u>SEE TABLE</u>			
<u>Drawn By</u>	<u>Approved By</u>			<u>Eyebolts</u>	<u>Mould Dimensions</u>			
Cassani M.	Direzione C.A.							
<u>N°Par.</u>	<u>N°Pez.</u>	<u>Dimensions</u>	<u>Material</u>	<u>Treatment</u>	<u>Scale</u>	A4		
						Sheet 1/1		

Changes
Signature
Date