



Luer cone ISO80369-7




I-I




Luer thread

Headquarters for ring



0.2 BURRS ALLOWED



MATERIAL	CODE
PC	CONREVCS
PC M	CONREVCS_M
ACRYLIC	CONREVCSY

 CARRO LUIGI MOLDS AND PLASTIC COMPONENTS Via Lavoro e industria 25030 Rudiano (BS) M0803.04 rev.01 del 02/02/2024 General tolerances ISO 2768-mk 0.5 < x ≤ 6 ± 0.1 6 < x ≤ 30 ± 0.2 30 < x ≤ 120 ± 0.3 120 < x ≤ 400 ± 0.5	Mould Denomination _____				N°Cavities ---	Order ---
	Member Denomination LULLABY REVOLVING HANDLE				Costumer code ---	
	Date 02/04/2025	Molding Material SEE TABLE	C.tion -----	Machine -----	Costumer part ---	
	Drawn By Lanza A.		Approved By Direzione C.A.		Costumer part ---	
	N°Par. ---	N°Pez. ---	Dimensions -----	Material -----	Treatment -----	Scale 2:1
					A4	Sheet 1/1

3

2

20/05/2026 Coccaglio 1 Updated title block

Date Signature Changes